

# GDI2 antibody

Catalog Number: 110927

## Product name

GDI2 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

GDI2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Positive IHC detected in human pancreas cancer. Observed molecular weight by Western-blot: 48kd

## Preparation

This antibody was obtained by immunization of GDI2 recombinant protein (Accession Number: NM\_001494). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:10-1:100


## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110927(GDI2 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:110927(GDI2 antibody) at dilution of 1:10 (under 10x lens)



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:110927(GDI2 antibody) at dilution of 1:10 (under 40x lens)