

H5N1 Neuraminidase (Active) recombinant protein



Catalog Number: 501927

General Information

Protein Construction

A DNA sequence encoding the Influenza A virus (A/Egypt/2321-NAMRU3/2007(H5N1)) neuraminidase (ACJ53852.1) (Met1-Lys449) was expressed, the cell lysates are collected, and bio-activity was tested.

Organism

H5N1

Expression Host

Human Cells

QC Testing

Activity

Measured by its ability to cleave a fluorogenic substrate, 2'-(4-Methylumbelliferyl)- α -D-N-acetylneuraminic acid The specific activity is > 100 U

One unit is defined as the amount of enzyme required to cleave 1 nmole of 2'-(4-Methylumbelliferyl)- α -D-N-acetylneuraminic acid per minute at pH 7.5 at 37°C.

Endotoxin

< 1.0 EU per μ g of the protein as determined by

the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Molecular Mass

The recombinant influenza A H5N1 Neuraminidase comprises 449 amino acids and has a predicted molecular mass of 49.1 kDa.

Formulation

Lyophilized from sterile PBS, 1% Triton X-100, 5% Trehalose, 5% Mannitol, pH7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide

Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.