

TMEM43 Polyclonal Antibody



Catalog Number: 165927

Product name

TMEM43 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to TMEM43

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 80-310 of human TMEM43 (NP_077310.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:100

Validations



Western blot - TMEM43 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using TMEM43 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 5s.



Immunohistochemistry - TMEM43 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using TMEM43 antibody at dilution of 1:100 (40x lens).