

# Anti-ALKBH2 Rabbit antibody



Catalog Number: 166927

## Product name

Anti-ALKBH2 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

ALKBH2 Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human ALKBH2

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, FC

## Dilutions

WB: 1/1000

FC: 1/20-1/100

## Validations



Western blot detection of ALKBH2 in C6,3T3,HeLa cell lysates using ALKBH2 Rabbit pAb(1:1000 diluted).Predicted band size:29kDa.Observed band size:29kDa.