

ZC3H7B antibody



Catalog Number: 116927

Product name

ZC3H7B antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ZC3H7B Rabbit Polyclonal antibody. Positive WB detected in HUVEC cells, mouse brain tissue. Observed molecular weight by Western-blot: 111 kDa

Preparation

This antibody was obtained by immunization of ZC3H7B recombinant protein (Accession Number: NM_017590). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HUVEC cells were subjected to SDS PAGE followed by western blot with Catalog No:116927(ZC3H7B Antibody) at dilution of 1:600