

DHX15 antibody

Catalog Number: 109927

Product name

DHX15 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DHX15 Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells, HEK-293 cells, HeLa cells. Positive IP detected in NIH/3T3 cells. Observed molecular weight by Western-blot: 91 kDa

Preparation

This antibody was obtained by immunization of DHX15 recombinant protein (Accession Number: NM_001358). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:109927(DHX15 antibody) at dilution of 1:300



IP Result of anti-DHX15 (IP:Catalog No:109927, 4ug; Detection:Catalog No:109927 1:600) with NIH/3T3 cells lysate 2000ug.