

EHHADH antibody



Catalog Number: 110227

Product name

EHHADH antibody

Immunogen

[Human EHHADH Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

EHHADH Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, mouse liver tissue, SMMC-7721 cells. Observed molecular weight by Western-blot: 69 kDa

Preparation

This antibody was obtained by immunization of EHHADH recombinant protein (Accession Number: NM_001966). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:1000-1:10000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:110227(EHHADH Antibody) at dilution of 1:2000