

Human PLTP (His Tag) recombinant protein



Catalog Number: 502227

General Information

Gene Name Synonym

Lipid transfer protein II

Protein Construction

A DNA sequence encoding the human PLTP isoform 1 (P55058-1) (Met 1-Val 493) was expressed, fused with a polyhistidine tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 80 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Glu 18

Molecular Mass

The secreted recombinant human PLTP consists of

487 amino acids after removal of the signal peptide and has a predicted molecular mass of 54.5 kDa. As a result of glycosylation, the apparent molecular mass of rh PLTP is approximately 70-80 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human PLTP Protein (His Tag) SDS-PAGE