

# TRA2B antibody



Catalog Number: 116227

## Product name

TRA2B antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

TRA2B Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, HepG2 cells. Observed molecular weight by Western-blot: 34 kDa

## Preparation

This antibody was obtained by immunization of TRA2B recombinant protein (Accession Number: NM\_004593). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:116227(TRA2B Antibody) at dilution of 1:1000