

Anti-CDKN2A/p16INK4a Rabbit antibody



Catalog Number: 167227

Product name

Anti-CDKN2A/p16INK4a Rabbit antibody

Specificity

Human

Antibody description

CDKN2A/p16INK4a Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human CDKN2A/p16INK4a

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/1000


IHC: 1/20

ICC/IF: 1/50

FC: 1/50

IP: 1/20

Validations



Western blot detection of p16 INK4A in 293 cell lysates using p16 INK4A Rabbit pAb(1:1000 diluted). Predicted band size: 17KDa. Observed band size: 17KDa.