

Human SEMA4B/KIAA1745/SEMAC (His tag) recombinant protein



Catalog Number: 518227

General Information

Gene Name Synonym

Semaphorin-4B, KIAA1745, SEMAC and SEMA4B

Protein Construction

A DNA sequence encoding the Human SEMA4B/KIAA1745/SEMAC Leu39-Glu712 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM PB, 150mM NaCl, pH7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human SEMA4B/KIAA1745/SEMAC (His tag)
recombinant protein