

IFNAR2 antibody



Catalog Number: 111627

Product name

IFNAR2 antibody

Specificity

Human; other species not tested.

Antibody description

IFNAR2 Rabbit Polyclonal antibody. Positive FC detected in Raji cells. Positive WB detected in MCF7 cells. Observed molecular weight by Western-blot: 58-66 kDa

Preparation

This antibody was obtained by immunization of IFNAR2 recombinant protein (Accession Number: NM_000874). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, FC

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:111627(IFNAR2 antibody) at dilution of 1:500



1X10⁶ Raji cells were stained with 0.2ug IFNAR2 antibody (Catalog No:111627, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.