

Human TFAP2C / AP2-GAMMA (His Tag) recombinant protein



Catalog Number: 502627

General Information

Gene Name Synonym

Activating enhancer-binding protein 2 gamma;
Transcription factor ERF-1

Protein Construction

A DNA sequence encoding the mature form of human TFAP2C (Q92754) (Leu128-Val223) was expressed with a polyhistidine tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

> 98 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant human TFAP2C consists of 111 amino acids and predicts a molecular mass of 12.3 KDa. It migrates as an approximately 12-14 KDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human TFAP2C / AP2-GAMMA Protein (His Tag)
SDS-PAGE