

PCGF6 Polyclonal Antibody



Catalog Number: 164627

Product name

PCGF6 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to PCGF6

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 161-350 of human PCGF6 (NP_001011663.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:100

Validations



Western blot - PCGF6 Polyclonal Antibody

Western blot analysis of extracts of MCF7 cells, using PCGF6 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.



Immunofluorescence - PCGF6 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using PCGF6 antibody.