

# ST5 antibody

Catalog Number: 115627

---

## Product name

ST5 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

ST5 Rabbit Polyclonal antibody. Positive WB detected in human testis tissue, HeLa cells, human adrenal gland tissue, mouse uterus tissue, SH-SY5Y cells. Observed molecular weight by Western-blot: 82kd,70kd

## Preparation

This antibody was obtained by immunization of ST5 recombinant protein (Accession Number: NM\_001376496). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:115627(ST5 antibody) at dilution of 1:500