

# TMEM106B antibody



Catalog Number: 107627

## Product name

TMEM106B antibody

## Specificity

Human; other species not tested.

## Antibody description

TMEM106B Mouse Monoclonal antibody. Positive WB detected in U-251 cells. Observed molecular weight by Western-blot: ~40 kDa

## Preparation

This antibody was obtained by immunization of TMEM106B recombinant protein (Accession Number: NM\_001134232). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



U-251 cells were subjected to SDS PAGE followed by western blot with Catalog No:107627(TMEM106B Antibody) at dilution of 1:1000