

ARMS2 Polyclonal Antibody



Catalog Number: 160327

Product name

ARMS2 Polyclonal Antibody

Specificity

Mouse

Antibody description

Polyclonal antibody to ARMS2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-107 of human ARMS2 (NP_001093137.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:200 - 1:1000

Validations



Western blot - ARMS2 Polyclonal Antibody

Western blot analysis of extracts of mouse eye, using ARMS2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.