

# Anti-SAA4 antibody



Catalog Number: 142327

## Product name

Anti-SAA4 antibody

## Immunogen

[Human SAA4 \(GST Tag\) recombinant protein](#)

## Specificity

Human SAA4

## Antibody description

Mouse Monoclonal to Human SAA4

## Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human SAA4 and conjugated with PE under optimum conditions, the unreacted PE was removed.

## Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

## Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from

prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

FCM

## Dilutions

FCM

## Validations



SAA4 Antibody, Mouse MAb, Flow Cytometry