

Anti-GFER / ALR antibody



Catalog Number: 106327

Product name

Anti-GFER / ALR antibody

Immunogen

[Human GFER / ALR \(His Tag\) recombinant protein](#)

Specificity

Human GFER / ALR

Antibody description

Rabbit polyclonal to GFER / ALR

Preparation

Produced in rabbits immunized with purified, recombinant Human GFER / ALR (rh GFER / ALR; NP_005253.3; Met1-Asp125). GFER / ALR specific IgG was purified by Human GFER / ALR affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA: 0.1-0.2 µg/ml

This antibody can be used at 0.1-0.2 µg/ml with the appropriate secondary reagents to detect Human GFER / ALR.

IHC-P: 0.1-2 µg/mL

Validations



GFER / ALR Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human GFER in human liver with rabbit polyclonal antibody (1 µg/mL, formalin-fixed paraffin embedded sections).