

Anti-Cardiac Troponin I (5C6) Mouse antibody



Catalog Number: 168327

Product name

Mouse IgG2b

Anti-Cardiac Troponin I (5C6) Mouse antibody

Applications

Specificity

WB, ELISA

Human, Mouse, Rat

Dilutions

Antibody description

ELISA: 1/10000

Cardiac Troponin I (5C6) Mouse monoclonal antibody

;WB: 1/1000

Preparation

Validations

Antigen: Synthetic peptides corresponding to the sequence of human Cardiac Troponin I.

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300.

Storage

Store at 4°C short term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

