

Anti-Ube2N/Ubc13 Rabbit antibody



Catalog Number: 167727

Product name

Anti-Ube2N/Ubc13 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Ube2N/Ubc13 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Ube2N/Ubc13

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF


Dilutions

WB: 1/1000

IHC: 1/50

ICC/IF: 1/50

Validations



Western blot detection of Ube2N/Ubc13 in HL-60, C6 using Ube2N/Ubc13 Rabbit mAb(1:1000 diluted)