

C4 (alpha chain) antibody



Catalog Number: 108727

Product name

C4 (alpha chain) antibody

Specificity

Human; other species not tested.

Antibody description

C4 (alpha chain) Rabbit Polyclonal antibody.
Positive WB detected in human blood tissue.
Observed molecular weight by Western-blot: 95 kDa

Preparation

This antibody was obtained by immunization of C4 (alpha chain) recombinant protein (Accession Number: NM_007293). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



human blood tissue were subjected to SDS PAGE followed by western blot with Catalog No:108727(C4A Antibody) at dilution of 1:3000