

GLT25D1 antibody



Catalog Number: 111037

Product name

GLT25D1 antibody

IP: 1:200-1:2000

IF: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GLT25D1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, SH-SY5Y cells. Positive IP detected in HeLa cells. Positive IF detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 72 kDa

Preparation

This antibody was obtained by immunization of GLT25D1 recombinant protein (Accession Number: BC020492). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111037(GLT25D1 Antibody) at dilution of 1:1000



IP Result of anti-GLT25D1 (IP:Catalog No:111037, 4ug; Detection:Catalog No:111037 1:500) with HeLa cells lysate 4000ug.



Immunofluorescent analysis of SH-SY5Y cells using Catalog No:111037(GLT25D1 Antibody) at

GLT25D1 antibody



Catalog Number: 111037

dilution of 1:50 and Alexa Fluor 488-conjugated

AffiniPure Goat Anti-Rabbit IgG(H+L)