

SAA1 antibody

Catalog Number: 115037

Product name

SAA1 antibody

Specificity

Human; other species not tested.

Antibody description

SAA1 Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue.

Preparation

This antibody was obtained by immunization of SAA1 recombinant protein (Accession Number: NM_199161). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC

Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations

Immunohistochemical of paraffin-embedded human liver using Catalog No:115037(SAA1 antibody) at dilution of 1:50 (under 40x lens)