

THAP1 antibody

Catalog Number: 116037

Product name

THAP1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

THAP1 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive WB detected in HEK-293 cells, A375 cells, HeLa cells, HepG2 cells, human brain tissue, human heart tissue, mouse brain tissue, mouse heart tissue, SGC-7901 cells. Observed molecular weight by Western-blot: 30kd

Preparation

This antibody was obtained by immunization of THAP1 recombinant protein (Accession Number: NM_018105). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:116037(THAP1 antibody) at dilution of 1:50



(under 40x lens)

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116037(THAP1 Antibody) at dilution of 1:1000