

H9N5 HA (His Tag) recombinant protein



Catalog Number: 501437

General Information

Protein Construction

A DNA sequence encoding the Influenza A virus (A/shorebird/DE/261/2003(H9N5)) hemagglutinin (ABB87950.1) (Met1-Lys523), termed as HA, was expressed with a polyhistidine tag at the C-terminus.

Organism

H9N5

Expression Host

Baculovirus-Insect Cells

QC Testing

Activity

Measured by its ability to agglutinate guinea pig red blood cells.

HA titer is 0.02-0.1 μ g/mL for 1%GRBC

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μ g protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 19

Molecular Mass

The recombinant hemagglutinin of Influenza A virus (A/shorebird/DE/261/2003(H9N5)) consists 516 amino acids and predicts a molecular mass of 58.3 kDa.

Formulation

Lyophilized from sterile 20 mM Tris, 500 mM NaCl, 10 % glycerol, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H9N5 (A/shorebird/DE/261/2003)
Hemagglutinin / HA Protein (His Tag)