

MAPK7 antibody



Catalog Number: 112437

Product name

MAPK7 antibody

Specificity

human, mouse; other species not tested.

Antibody description

MAPK7 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, HeLa cells, human heart tissue, mouse brain tissue. Observed molecular weight by Western-blot: 59kd

Preparation

This antibody was obtained by immunization of MAPK7 recombinant protein (Accession Number: XM_011523957). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:112437(MAPK7 antibody) at dilution of 1:300