

Anti-IGSF8 antibody



Catalog Number: 176437

Product name

Polyclonal

Anti-IGSF8 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse

Applications

Antibody description

WB

Rabbit polyclonal antibody to IGSF8

Dilutions

Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human igsf8. aa range:451-500

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Validations




Fig1.; Western Blot analysis of NIH-3T3 cells using CD316 Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000