

# Anti-Prostaglandin D Synthase Rabbit antibody



Catalog Number: 167437

## Product name

Rabbit IgG

Anti-Prostaglandin D Synthase Rabbit antibody

## Applications

## Specificity

WB, IHC-P, ICC/IF

Human, Mouse, Rat

## Dilutions

## Antibody description

WB: 1/1000

Prostaglandin D Synthase Rabbit polyclonal antibody

IHC: 1/20

## Preparation


ICC/IF: 1/20

Antigen: A synthetic peptide of human Prostaglandin D Synthase

## Validations

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

Western blot detection of Prostaglandin D Synthase in K562, C6, 3T3, HeLa cell lysates using Prostaglandin D Synthase Rabbit pAb(1:1000 diluted). Predicted band size: 21kDa. Observed band size: 26kDa.

## Ig Type