

Anti-SOD2 (4C4) Mouse antibody



Catalog Number: 168437

Product name

Anti-SOD2 (4C4) Mouse antibody

Ig Type

Mouse IgG2b

Specificity

Human

Applications

WB

Antibody description

SOD2 (4C4) Mouse monoclonal antibody

Dilutions

WB: 1/1000

Preparation

Antigen: Synthetic peptide corresponding to human SOD2 protein

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Validations



Western blot detection of SOD2 in HeLa, THP-1, A549, Ramos and K562 cell lysates using SOD2 Mouse mAb (1:1000 diluted). Predicted band size: 25 kDa. Observed band size: 22 kDa.