

Collagen Type IV antibody



Catalog Number: 109437

Product name

Collagen Type IV antibody

Specificity

Human; other species not tested.

Antibody description

Collagen Type IV Rabbit Polyclonal antibody.
Positive IHC detected in human kidney tissue.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001846).
Purification method: Antigen affinity purified.

Formulation

0.1M NaHCO₃, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:109437(COL4A2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:109437(COL4A2 Antibody) at dilution of 1:50 (under 40x lens)