

# MGAT5B Polyclonal Antibody



Catalog Number: 161837

## Product name

MGAT5B Polyclonal Antibody

## Specificity

Mouse

## Antibody description

Polyclonal antibody to MGAT5B

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 616-790 of human MGAT5B (NP\_653278.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



Western blot - MGAT5B Polyclonal Antibody

Western blot analysis of extracts of mouse brain, using MGAT5B antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.