

Human USP14/TGT (His tag) recombinant protein



Catalog Number: 518837

General Information

Gene Name Synonym

Ubiquitin Carboxyl-Terminal Hydrolase 14,
Deubiquitinating Enzyme 14, Ubiquitin
Thioesterase 14, Ubiquitin-Specific-Processing
Protease 14, USP14, TGT

Protein Construction

A DNA sequence encoding the Human USP14/TGT
Asp91-Gln494 is expressed with a 6His tag at the
N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL
method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM
TrisHCl, 100mM NaCl, 20% Glycerol, pH 8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human USP14/TGT (His tag) recombinant protein