

# DDX51 antibody



Catalog Number: 109837

## Product name

DDX51 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

DDX51 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, Jurkat cells, Raji cells. Observed molecular weight by Western-blot: 72kd

## Preparation

This antibody was obtained by immunization of DDX51 recombinant protein (Accession Number: NM\_175066). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:109837(DDX51 antibody) at dilution of 1:300