

# DYNC1LI1 antibody



Catalog Number: 110137

## Product name

DYNC1LI1 antibody

## Immunogen

[Human DYNC1LI1 Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

DYNC1LI1 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive WB detected in K-562 cells, HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 57-60 kDa

## Preparation

This antibody was obtained by immunization of DYNC1LI1 recombinant protein (Accession Number: NM\_016141). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:110137(DYNC1LI1 Antibody) at dilution of



1:1000 Immunofluorescent analysis of (-20oc Ethanol) fixed HeLa cells using Catalog No:110137(DYNC1LI1 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)