

MYCT1 Polyclonal Antibody



Catalog Number: 162137

Product name

MYCT1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to MYCT1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 99-235 of human MYCT1 (NP_079383.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

Validations




Western blot - MYCT1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using MYCT1 antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunofluorescence - MYCT1 Polyclonal Antibody

Immunofluorescence analysis of C6 cells using MYCT1 antibody at dilution of 1:100. Blue: DAPI



for nuclear staining.

MYCT1 Polyclonal Antibody



Catalog Number: 162137

Immunofluorescence - MYCT1 Polyclonal Antibody
Immunofluorescence analysis of HeLa cells using

MYCT1 antibody at dilution of 1:100. Blue: DAPI
for nuclear staining.