

NF2 antibody



Catalog Number: 113137

Product name

NF2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NF2 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, human kidney tissue, Jurkat cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 64-66 kDa

Preparation

This antibody was obtained by immunization of NF2 recombinant protein (Accession Number: BC007279). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:113137(NF2 Antibody) at dilution of 1:600