

Anti-KAP1 Rabbit antibody



Catalog Number: 167137

Product name

Anti-KAP1 Rabbit antibody

Specificity

Human, Monkey

Antibody description

KAP1 Rabbit polyclonal antibody

Preparation

Antigen: Peptide sequence around aa.827~831 (L-S-G-G-P) derived from Human TIF1 β .

Formulation

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB:1/1000

Validations



Western blot analysis of extract from 3T3,PC12,COS7 and HeLa cells using TIF1u03b2 rabbit pAb (dilution 1:1000).Predicted band size:110KDa.Observed band size:110KDa.