

# Anti-FSH beta Mouse antibody



Catalog Number: 169137

## Product name

Anti-FSH beta Mouse antibody

## Specificity

Human

## Antibody description

FSH beta Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of human FSHB (AA: 19-129) expressed in E. Coli.

## Formulation

Purified antibody in PBS with 0.05% sodium azide.

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal


## Ig Type

Mouse IgG1


## Applications

WB, FC

## Validations




Western blot analysis using FSHB mAb against



human FSHB (AA

Flow cytometric analysis of HepG2 cells using FSHB mouse mAb (green) and negative control



(red).


Western blot analysis using FSHB mAb against

# Anti-FSH beta Mouse antibody



Catalog Number: 169137

---



Western blot analysis using FSHB mouse mAb against Cos7 (1), HepG2 (2), HEK293 (3) cell lysate.

HEK293 (1) and FSHB (AA