

P15RS antibody

Catalog Number: 113537

Product name

P15RS antibody

Immunogen

[Human P15RS Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

P15RS Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, A431 cells, Jurkat cells, L02 cells, mouse liver tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 36 kDa

Preparation

This antibody was obtained by immunization of P15RS recombinant protein (Accession Number: NM_001303413). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:113537(P15RS antibody) at dilution of 1:1000



Immunofluorescent analysis of HepG2 cells, using RPRD1A antibody Catalog No:113537 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).