

# BTBD6 antibody

Catalog Number: 108537

## Product name

BTBD6 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

BTBD6 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells. Positive IP detected in Y79 cells. Observed molecular weight by Western-blot: 53-58 kDa

## Preparation

This antibody was obtained by immunization of BTBD6 recombinant protein (Accession Number: BC042525). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IP, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:108537(BTBD6 Antibody) at dilution of 1:300



IP Result of anti-BTBD6 (IP:Catalog No:108537, 4ug; Detection:Catalog No:108537 1:300) with Y79 cells lysate 800ug.