

Anti-CYP21A2 antibody



Catalog Number: 175937

Product name

Anti-CYP21A2 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to CYP21A2

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human cytochrome p450 21a2. aa range:151-200

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IF

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

IF: 1:100-1:500

Validations




Fig1::; Western Blot analysis of 293 cells using CYP21A2 Polyclonal Antibody diluted at 1: 1000