

LIN7B antibody



Catalog Number: 112237

Product name

LIN7B antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

LIN7B Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, A549 cells, human brain tissue, K-562 cells. Observed molecular weight by Western-blot: 23kd

Preparation

This antibody was obtained by immunization of LIN7B recombinant protein (Accession Number: NM_022165). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:112237(LIN7B antibody) at dilution of 1:500