

NKX3-2 antibody



Catalog Number: 113237

Product name

NKX3-2 antibody

Immunogen

[Human NKX3-2 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

NKX3-2 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, mouse small intestine tissue. Observed molecular weight by Western-blot: 36-40 kDa

Preparation

This antibody was obtained by immunization of NKX3-2 recombinant protein (Accession Number: NM_001189). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:113237(NKX3-2 Antibody) at dilution of 1:300