

# Anti-OR9Q1 antibody

Catalog Number: 176237

---

## Product name

Anti-OR9Q1 antibody

## Specificity

Human

## Antibody description

Rabbit polyclonal antibody to OR9Q1

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or9q1. aa range:234-283

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IF

## Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

## Validations




Fig1:; Western Blot analysis of Jurkat cells using Olfactory receptor 9Q1 Polyclonal Antibody