

## General Information

### Gene Name Synonym

Alpha-2-antiplasmin, Alpha-2-AP, Alpha-2-plasmin inhibitor, Alpha-2-PI, Serpin F2, Serpinf2

### Protein Construction

A DNA sequence encoding the Mouse Serpinf2/Pli Val28-Lys491 is expressed with a Fc tag at the C-terminus.

### Organism

Mouse

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL

method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation

Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM Tris, 150mM NaCl, pH 8.0 (20% glycerol, 1mM DTT).


## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE

kDa MK R NR



Mouse Serpinf2/Pli (Fc tag) recombinant protein