

# PDS5A antibody



Catalog Number: 113637

## Product name

PDS5A antibody

IHC: 1:20-1:200

IF: 1:10-1:100


## Specificity

Human, Mouse, Rat; other species not tested.

## Validations

## Antibody description

PDS5A Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IHC detected in human colon tissue. Positive WB detected in HeLa cells, COLO 320 cells, HEK-293 cells, mouse kidney tissue. Observed molecular weight by Western-blot: 130kd




## Preparation

This antibody was obtained by immunization of PDS5A recombinant protein (Accession Number: BC009650). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113637(PDS5A antibody) at dilution of 1:400

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal


Immunofluorescent analysis of HeLa cells, using PDS5A antibody Catalog No:113637 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF, IHC



## Dilutions

Recommended Dilution:

WB: 1:500-1:5000


Immunohistochemistry of paraffin-embedded

# PDS5A antibody

Catalog Number: 113637

---

human colon tissue slide using Catalog No:113637(PDS5A Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human colon tissue slide using Catalog No:113637(PDS5A Antibody) at dilution of 1:50 (under 40x lens)