

# RGMA antibody

Catalog Number: 114637

## Product name

RGMA antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RGMA Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Positive WB detected in mouse testis tissue. Observed molecular weight by Western-blot: 49 kDa

## Preparation

This antibody was obtained by immunization of RGMA recombinant protein (Accession Number: BC015886). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:114637(RGMA antibody) at dilution of 1:400



IP Result of anti-RGMA (IP:Catalog No:114637, 4ug; Detection:Catalog No:114637 1:400) with mouse testis tissue lysate 4000ug.