

STAC2 antibody

Catalog Number: 115637

Product name

STAC2 antibody

Specificity

Human; other species not tested.

Antibody description

STAC2 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells. Positive IHC detected in human prostate hyperplasia tissue. Observed molecular weight by Western-blot: 42-45 kDa

Preparation

This antibody was obtained by immunization of STAC2 recombinant protein (Accession Number: NM_001351360). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:115637(STAC2 Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human prostate hyperplasia slide using Catalog No:115637(STAC2 Antibody) at dilution of 1:50