

# Anti-GC1q R Rabbit antibody



Catalog Number: 167637

## Product name

Anti-GC1q R Rabbit antibody

## Specificity

Human

## Antibody description

GC1q R Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human GC1q R

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, FC

## Dilutions


WB: 1/1000

IHC: 1/50

ICC/IF: 1/50

FC: 1/50

## Validations



Western blot detection of GC1q R in A549,HL-60 cell lysates using GC1q R Rabbit mAb(1:1000 diluted).Predicted band size:31kDa.Observed band size:31kDa.